



### STEP 1

Insert one Floss Capsule into the nozzle



### STEP 4

Press On/Off button to let the water blow the floss out into nozzle (prime the floss)



### STEP 2

Place the nozzle into the handle completely



### STEP 5

Put the nozzle toward your teeth and slightly close your lips. Lean over the sink and press On/Off button to start cleaning



### STEP 3

Connect the water supply hose to the handle and fill cup with water



Use JETPIK UV Sanitizer regularly to kill bacteria and keep nozzle and toothbrush clean

JETPIK web site: [www.jetpik.com](http://www.jetpik.com)

## Quick Start Guide

## JETPIK Power Floss System



JETPIK web site: [www.jetpik.com](http://www.jetpik.com)